

LED CELL ACTIVATOR (NB-306LCA)

LED CELL ACTIVATOR is effective for activating cell or growth factor before using.

Features

- Selective wavelength according to the purpose of treatment.
- RED : to help cell's growth, collagen synthesis
- BLUE : to help cell's self-purification, to restrain lipolysis
- GREEN : to help cellular immunity increase.
- YELLOW : to stimulate cell's activation.
- N-BIOTEK'S Patent product.



Specification



Items	Unit	NB-306LCA
LED Light		Red/Yellow/Green/Blue
Time		
range		Continuous or up to 47h 59min.
accuracy		± 1%
increment		1 Minute
Control		Microprocessor digital PID
Display		LED Display
Operating Panel		Touch Button
Block Material		Solid anodized aluminum
Block Capacity		50ml Conical Tube 1Holex2ea
Dimension (in)		30 Ø x100mmx2ea
(out)	mm	250(W)x250(D)x175(H)mm
Power	V/Hz	110/220V, 50/60Hz, 30W
Weight	kg	5kg